



# Jason<sup>®</sup> membrane

**NATIVE PERICARDIUM** GBR/GTR MEMBRANE

barrier membrane



native

reliable

resorbable

# botiss regeneration system



## Development / Production / Distribution

maxresorb®	maxresorb® inject	cerabone®	cerabone® plus	maxgraft®	maxgraft® cortico	maxgraft® bonering	maxgraft® bonebuilder
Synthetic biphasic calcium phosphate	Synthetic injectable bone paste	100% pure bovine bone mineral	cerabone® mixed with hyaluronate	Processed allogenic bone graft	Processed allogenic bone plate	Processed allogenic bone ring	Patient matched allogenic bone implant
NOVAMag® fixation screw	NOVAMag® membrane	permamem®	Jason® membrane	collprotect® membrane	mucoderm®	collacone®	collafleece®
Resorbable magnesium screw	Resorbable magnesium membrane	High-density PTFE barrier membrane	Native pericardium GBR / GTR membrane	Native collagen membrane	3D-stable soft tissue graft (Collagen)	Collagen hemostat (Cone)	Collagen hemostat (Sponge)

## 360° – the botiss regeneration system: Innovation, Safety, Reliability, and Aesthetics

**botiss biomaterials offers a unique systematic BTR approach - the complete regenerative biomaterial portfolio for Implantology, Oral and CMF Surgery, and Periodontology at hand.**

We all know – no single bone graft or soft tissue biomaterial can suit all medical needs, biological situations, and indications. Factors, such as indication, age, hygiene, biotype, bone height, and treatment plan, require a sophisticated approach with different, coordinated products.

To achieve optimal results, we offer you the botiss regeneration system. It includes all long-term proven biological materials (e.g., bovine, synthetic, allografts, collagen, granules, blocks, membranes, and soft tissue matrices), which can be used in various combinations for each specific indication. All products are manufactured according to the highest quality standards.

Patient's safety, ease of use and reliable treatment results – these are your and our first priorities. The products of the botiss regeneration system have proven their success in terms of safety, efficacy, and reliability in a multitude of preclinical and clinical studies and, most importantly, in the daily clinical work, with hundreds of thousands of patients treated worldwide.

We substantially invest in research and education. Unique innovations, such as cerabone® plus and NOVAMag®, the concept of high-quality learning and education with the botiss academy, and our international bone & tissue

days are the results of our partnership with world-wide renowned academic research institutes, global opinion leaders, and practitioners in their daily clinical environment.

botiss biomaterials is one of the leading companies in the field of dental bone and tissue regeneration. The botiss regeneration system is available in over 100 countries worldwide via a global network of distribution partners and employees, who are committed experts in the field of oral surgery and implantology.

botiss biomaterials is an innovative, clinically oriented medical device/pharmaceutical company headquartered in Germany and further development and production sites in Germany, Austria and England.

We proudly welcome you to the botiss regeneration system community. We invite you to share your experiences and suggestions with us, which are precious to further improve our products or develop new product concepts.

Dr. Drazen Tadic  
[dt@botiss.com](mailto:dt@botiss.com)


Oliver Bielenstein  
[ob@botiss.com](mailto:ob@botiss.com)



# Jason® membrane

NATIVE PERICARDIUM GBR/GTR MEMBRANE

Jason® membrane is a native collagen membrane obtained from porcine pericardium, developed and manufactured for dental bone and tissue regeneration. Besides its very low thickness it provides a naturally long barrier function based on the specific composition and structure of the pericardial collagen fibers<sup>7,9</sup>.



SEM image of Jason® membrane

Owing to the preservation of the natural biomechanical properties of the pericardium, Jason® membrane exhibits beneficial handling characteristics such as a remarkable tear resistance and effective surface adaptation<sup>8</sup>.



## INDICATIONS:

Implantology,  
Periodontology and  
Oral and CMF Surgery

- Fenestration and dehiscence defects
- Sinus lift
- Socket and ridge preservation
- Alveolar ridge augmentation and reconstruction
- Intraosseous defects (1 to 3 walls)
- Furcation defects (class I and II)
- Protection and covering of minor perforations of the Schneiderian membrane



Jason® membrane maintained the barrier function 56 days after subcutaneous implantation in rats

Due to the unique production process, the superior properties of the native pericardium are preserved during the extensive cleaning procedure that is applied for the production of Jason® membrane. Therefore, Jason® membrane shows a natural honeycomb-like, multilayered collagen structure with an increased content of collagen type III leading to a remarkable tear resistance and slow degradation of Jason® membrane<sup>9, 10, 11, 12</sup>. This ensures a natural long barrier function, making the Jason® membrane our recommended choice particularly for large augmentative procedures.

## PROPERTIES

- Naturally long barrier function
- Multi-directional strength and tear resistance
- No stickiness after hydration
- Excellent surface adaptation
- Easy handling, can be applied dry or wet
- Low thickness

## PRODUCT SPECIFICATIONS

Art.No.	Size	Content
681520	15 x 20 mm	1 membrane
682030	20 x 30 mm	1 membrane
683040	30 x 40 mm	1 membrane



# THE ORIGIN OF COLLAGEN MEMBRANES

The first collagen membranes available on the market were of bovine origin (Achilles tendon and pericardium). Nowadays, porcine membranes are widely used because their application excludes the risk of BSE transmission.

Moreover, porcine collagen exhibits a high homology to human collagen and therefore a very low antigenicity<sup>6</sup>.


Due to these reasons, today collagen membranes are exclusively produced from porcine collagen.

Collagen membranes may be derived from various tissues, ranging from dermis, to peritoneum and pericardium. Accordingly, these membranes differ in their handling and degradation properties, as well as their barrier function.

## THE ADVANTAGES OF COLLAGEN

Several factors make collagen an optimal biological material for the use as barrier membranes. One important characteristic is the excellent biocompatibility of collagen and its degradation products. Collagen is widely distributed throughout the body, making up approx. 60% of all proteins within the gingival connective tissue.

Due to their low antigenicity, animal collagens may be used in humans without causing tissue rejection.



3D structure of a collagen fleece


Collagens are resistant to any unspecific proteolytic degradation and are only degraded by specific enzymes called collagenases. Collagens are involved in the primary hemostatic reaction. Thus, collagen membranes contribute to a fast stabilization of the wound area. Another advantage of collagen is its chemotactic attraction of regenerative cells such as osteoblasts, gingival fibroblasts and periodontal ligament cells. Following dehiscence, the exposure of a collagen membrane leads to its quick proteolytic degradation. However, a secondary granulation without any inflammatory reaction may be observed<sup>4</sup>.

### BENEFITS

- Exceptional biocompatibility
- Support of hemostasis
- Low antigenicity
- Degradation by specific enzymes
- Protection and covering of minor perforations of the Schneiderian membrane

## COLLAGEN – A NATURAL HEMOSTATIC AGENT

Damage to the blood vessel wall leads to subendothelial collagen exposure. The collagen directly or indirectly interacts with the surface receptors on thrombocytes. The binding of collagen initiates a reaction cascade leading to transformation and aggregation of the thrombocytes. Additionally, the thrombocytes are cross-linked by fibrinogen. The resulting (white) thrombus initially stabilizes the wound<sup>5</sup>. Accordingly, collagen membranes support the formation of a blood coagulum and contribute to a rapid stabilization of the wound area. Due to their hemostatic effect, collagens are not only used as barrier membranes, but also as collagen sponges and cones for stabilization of biopsy harvesting sites or covering of minor oral wounds and extraction sockets.



## PROPERTIES OF PERICARDIUM MEMBRANES

Many collagen membranes have a limited barrier function due to their rapid enzymatic degradation. The stability and barrier function of collagen membranes are tightly linked to the properties of the native tissue from which they originate.

Jason<sup>®</sup> membrane is a native collagen membrane **derived from porcine pericardium.**

The pericardium is a double-layered tissue sac that surrounds and protects the heart, and allows volume changes by contraction and relaxation of the heart muscle. Imagine the human heart beats around 100.000 times a day, that's almost 40 million times per year! For every beat, the heart requires a similar effort as that needed to squeeze a tennis ball with one hand. The pericardial tissue consists of highly organized, naturally cross-linked fibers of collagen type I/III. During the production process, this particular collagen structure and composition is preserved.



Due to its structural characteristics Jason<sup>®</sup> membrane undergoes slow degradation and thus offers a prolonged barrier function<sup>7</sup>. Furthermore, Jason<sup>®</sup> membrane is distinguished by extraordinarily high tear resistance and excellent handling properties (e.g. good adaptation to surface contours, no sticking)<sup>8</sup>.

Despite its low thickness, Jason<sup>®</sup> membrane exhibits an excellent multidirectional tear resistance

# Production process




Jason® membrane

## SAFETY AND STABILITY

Jason® membrane consists of natural porcine collagen originating from animals destined for the food industry and certified according to EN ISO 22442.

### PERICARDIUM

Jason® membrane is a native membrane, the natural properties of the original tissue (pericardium) is preserved during the production process. The inherent architecture of the collagen structure provides superior handling properties, such as tear resistance, tensile strength, and adaptation to surface contours, in comparison to „non-native“ collagen membranes<sup>10</sup>.



Natural three-dimensional collagen network of Jason® membrane

The particular multi-stage cleaning process effectively removes all non-collagenic proteins and antigenic components. The resulting membranes exhibit a natural three-dimensional collagen structure mainly composed of collagen type I and a lower share of collagen type III.


Multi-stage cleaning process

Lyophilization

Packaging  
EO-sterilization

Sterile product

Jason® membrane




# Pre-clinical testing

## JASON® MEMBRANE SUPPORTS ATTACHMENT AND PROLIFERATION OF OSTEOBLAST-LIKE CELLS


Results of *in vitro* cell cultures. Dr. M. Herten, University of Münster and Prof. Dr. D. Rothamel, Mönchengladbach Hospital, University of Düsseldorf<sup>7</sup>

Incubation of the multi-layered Jason® membrane and a competitive bilayer membrane with osteoblast-like SaOs-2 cells showed a significantly higher cell proliferation on the Jason® membrane after seven days. The excellent cell attachment and proliferation on Jason® membrane highlights its suitability as scaffold for osteoblast guidance, which supports the bony regeneration of covered defects.




## DEGRADATION AND BARRIER FUNCTION OF JASON® MEMBRANE

Results from a degradation study in a mouse model<sup>9</sup>



The indication, the way of application as well as differences in the metabolism of each patient affect the degradation rate of a collagen membrane. To follow the progression of resorption of Jason® membrane at different time points, animal experiments were performed as these tests cannot be realized in humans. Subcutaneous implantation in mice demonstrated that Jason® membrane was intact after a period of 60 days. As the metabolic rate of mice is much higher compared to humans, the barrier function following implantation in humans will be significantly longer.



Pictures and data adapted from Barbeck et al. 2014<sup>9</sup> showing the degradation of Jason® membrane after subcutaneous implantation in mice. Membrane is still intact after 60 days.


# In vivo pre-clinical testing

## Jason® membrane – EXCELLENT BIOCOMPATIBILITY AND TISSUE INTEGRATION


Results from an animal model, Prof. Dr. Dr. D. Rothamel, Mönchengladbach Hospital, University of Düsseldorf<sup>7</sup>

Analysis of the tissue integration and morphological structure of Jason® membrane at four to 12 weeks after lateral augmentation in a dog model.

The membrane was integrated into the surrounding tissue without any inflammation. Significant degradation of the membrane started at week eight and proceeded until week 12. A bilayer membrane that was tested in the same model showed a comparably good tissue integration, but was almost completely degraded after eight weeks<sup>7</sup>.




Jason® membrane after four weeks healing time



The bilayer membrane after four weeks healing time

**4 weeks** healing time  
Both membranes showed good tissue integration without any inflammatory reaction, as demonstrated by Toluidine staining. Initial ingrowth of blood vessels improves nutrition of the graft and osseous regeneration.

**8 weeks** healing time  
The bilayer membrane was almost completely resorbed. Jason® membrane was still intact, serving as barrier against ingrowth of surrounding soft tissue.



The bilayer membrane after eight weeks healing time



Jason® membrane after eight weeks healing time



Jason® membrane after 12 weeks healing time

**12 weeks** healing time  
Jason® membrane was almost completely degraded and replaced by a periosteum rich in collagen fibers. The collagen of the membrane is partially visible as cloudy fibrous areas.

## CLINICAL CASE BY Dr. Sebastian Stavar, Houten, Netherlands

### RIDGE AUGMENTATION




Initial clinical situation with broken bridge abutment in regio 12, tooth 21 not worth preserving and tooth 11 lost by a front teeth trauma several years ago



Situation after atraumatic tooth extraction and suturing of wound margins



Clinical situation five weeks after extraction



Preparation of a mucoperiosteal flap - extensive bone deficit in horizontal and vertical dimension



Horizontal and vertical augmentation with cerabone® and autologous bone after placement of two implants



Coverage of the augmentation site with Jason® membrane



Tension-free wound closure



Clinical view two weeks post-operative



Complication free healing eleven weeks after augmentation



Exposure of implants and insertion of healing abutments



Shaping of the emergence profile using the temporary prosthesis



Final prosthetic restoration with implant-borne bridge in regio 12-21 and crown on tooth 22

**CLINICAL CASE BY**

Prof. Dr. Dr. Daniel Rothamel, Mönchengladbach Hospital,  
University of Düsseldorf, Germany


**SINUS LIFT WITH TWO-STAGE IMPLANT PLACEMENT**



Clinical situation before sinus lift



Clinical situation before sinus lift, occlusal view




Clinical situation following preparation of the mucoperiosteal flap



Preparation of a lateral sinus window



Placing of Jason® membrane in the sinus cavity



Jason® membrane serves as protection for the Schneiderian membrane



Filling the sinus cavity with cerabone®



cerabone® in the sinus cavity



Additional lateral augmentation with cerabone®



Covering of the augmentation area with Jason® membrane




Tension-free wound closure with single interrupted sutures



Excellent osseous integration of the cerabone® particles without soft tissue ingrowth at re-entry, six months post-operative



Stable insertion of two implants into sufficient bone matrix



Histological sections of biopsy taken at the time of implantation



Magnification of the histological sample demonstrates complete integration of cerabone® particles within the newly formed bone matrix



Post-operative X-ray

**CLINICAL CASE BY**

Prof. Dr. Dr. Daniel Rothamel, Mönchengladbach Hospital,  
University of Düsseldorf, Germany

**LATERAL AUGMENTATION**



Lateral defect in regio 24 at six months after extraction



Crestal view of defect



Surgical presentation of the bone defect



Thin buccal bone after implant installation



Dehiscence defect at palatal side



Lateral augmentation with cerabone® and autologous bone (mixture 1:1)



Further augmentation at the palatal side



Application of Jason® membrane



Soft tissue closure



Clinical situation after three months



Satisfactory bone formation and volume maintenance



Stable hard tissue conditions on both buccal and palatal side


**CLINICAL CASE BY**

Dr. Dr. Dr. Oliver Blume, Munich, Germany


**RIDGE AUGMENTATION IN THE MAXILLA**



Preoperative clinical situation - severe atrophy of the maxillary bone



Three dimensional reconstruction of the bone defect and planned maxgraft® bonebuilder blocks (blue)




Upper left maxilla - severe atrophic ridge




Fixation of maxgraft® bonebuilder and contouring with allogenic particulated material




Covering with Jason® membrane and one layer of PRF matrices



Tension-free and saliva-proof wound closure



Fixation of two more maxgraft® bonebuilder blocks on upper right maxillary ridge



X-ray six months post-operative



Clinical situation six months after augmentation



Implant placement



Temporary provision


**CLINICAL CASE BY**

Dr. Marius Steigmann, Neckargemünd, Germany

**CERABONE® FOR COVERAGE OF IMPLANT DEHISCENCE AND RIDGE AUGMENTATION**



Extraction of tooth 21 after endodontic treatment



Application of collacone® for stabilization of the blood clot



Buccal bone defect after eight weeks healing time




A periodontal probe demonstrates the vertical extension of the defect




Implant placed into the former extraction socket



Surface of the implant is covered with autologous bone



Coverage of the autologous bone with cerabone® (0.5 - 1.0 mm)



Covering of the bone substitute with Jason® membrane



Closure of the site using single sutures after periosteum slitting



Tension-free suturing maintains undisturbed healing



Abutment installation after implant uncovering, six months after implantation



Final prosthetic restoration with a full-ceramic crown



Radiographic control five years post-operative



**CONTOUR MAINTENANCE**

For augmentations in the aesthetic region cerabone® provides long-term dimensional stability and therefore a good bone bed to support an optimal contour of the soft tissue and sustained aesthetic result.

bone & tissue  
regeneration

botiss  
biomaterials

# Innovation. Regeneration. Aesthetics.

soft tissue

education

hard tissue

botiss biomaterials GmbH  
Hauptstr. 28  
15806 Zossen  
Germany

Tel.: +49 33769 / 88 41 985  
Fax: +49 33769 / 88 41 986

[www.botiss.com](http://www.botiss.com)  
facebook: [botissdental](https://www.facebook.com/botissdental)  
instagram: [botiss\\_dental](https://www.instagram.com/botiss_dental)

## LITERATURE

1. Dahlin et al. (1988). Healing of bone defects by guided tissue regeneration. *Plast Reconstr Surg.* 81(5).
2. Scantlebury (1994). 1982-1992: A decade technology development for guided tissue regeneration. *Periodontol*; 64:1129-1137.
3. Rothamel et al. (2005). Biodegradation of differently cross-linked collagen membranes: an experimental study in the rat. *Clin Oral Implants Res* 16:369-378.
4. Schwarz et al. (2006). Einsatz nativer und quervernetzter Kollagenmembranen für die gesteuerte Gewebe- und Knochenregeneration. *SCHWEIZ MONATSSCHR ZAHNMED* 116(11): 1112.
5. Nuytens et al. (2011). Platelet adhesion to collagen. *Thromb Res* 127 Suppl 2:S26-9.
6. Silvipriya et al. (2015). Collagen: Animal Sources and Biomedical Application. *Journal of Applied Pharmaceutical Science* 5, 123-127.
7. Rothamel et al. (2012). Biocompatibility and Biodegradation of a Native, Porcine Pericardium Membrane. Results from in vitro/in vivo Examination. *J. Int J Oral Maxillofac Implants* 2012, 27(1):146-54.
8. Usability testing Jason® membrane, data on file
9. Barbeck et al. (2015). Porcine Dermis and Pericardium-Based, Non-Cross-Linked Materials Induce Multinucleated Giant Cells After Their In Vivo Implantation: A Physiological Reaction? *J Oral Implantol.* 41(6):e267-81.
10. Braga-Vilela et al. (2008). Extracellular matrix of porcine pericardium: biochemistry and collagen architecture. *J Membrane Biol.* 221(1):15-25.
11. Pabst et al. (2020). Collagen matrices: opportunities and perspectives in oral hard and soft tissue regeneration. *Quintessence Int.*, 51(4):318-327.
12. Naimark et al. (1992). Correlation of structure and viscoelastic properties in the pericardia of four mammalian species. *Am J Physiol.* 1992 Oct;263(4 Pt 2):H1095-106.